



Colpotomy set

Uncomplicated use in gynecology

Reusable and universal instrument to be used in colpotomy, LAVH or TLH.

1 Atraumatic instrument attachment

In gynecologic-anatomy adapted design for use with ultracision

2 Optimal sizes

Reusable caps in three sizes

5 Stable access

One-way valve as an access for 10 mm instruments

3 Simple reprocessing

robust, biocompatible high performance plastic can be easily processed at 134 ° C / 3 bar

4 Lightweight

solid, reusable basic tube



- Good visualization of the uterine vessels for a safe dissection
- Stable pneumoperitoneum during application
- Ease of use, increases patient safety and allows a reduced operating time