



Light Solution

high performance plastic combined with high-quality stainless steel

Light and Reusable

The LightSolution product range is made of high performance plastic housing combined with high quality stainless steel tube, optimal working length and unique colour coding.


The solution for minimally invasive procedures

1 Different working lengths

- from 110 to 120 mm

2 Reducer included

- from Ø 11 mm
down to Ø 5,5 mm
- from Ø 13 mm
down to Ø 5,5 mm



3 Finger flange

- for better handling

4 Tube surfaces

- threaded tube

5 Colour coding

- ○ Ø 5,5 mm
- ● Ø 11 mm
- ● Ø 13 mm

- Products available in different diameters and working lengths
- Finger flange for better handling