



Light Solution high-performance plastic

Light and Reusable

The LightSolution product range is made of high-performance plastic, with optimal working length and unique colour coding.


The solution for minimally invasive procedures

1 Different working lengths

- from 120 to 130 mm

2 Reducer included

- from Ø 11 mm
down to Ø 5,5 mm
- from Ø 13 mm
down to Ø 5,5 mm



3 Finger flange

- for better handling

4 Tube surfaces

- threaded tube

5 Colour coding

- ○ Ø 5,5 mm
- ● Ø 11 mm
- ● Ø 13 mm

- Products available in different diameters and working lengths
- Finger flange for better handling