

Convincing quality in every detail

Top form in performance and function!

Controlled force – high torque


The special full load allows a controlled, powerful cutting of tissue with a very low drive.

Safe application –

The activating ring with integrated „Cut / No-Cut“ function increases the safety and protects the blade during nonapplication.

Pressure resistant –

The special Valve Module enables the use of Trokamed Tenaculum Forceps and maintains pressure constantly even during the change of instruments.



Sharp cut with big lumen –

The Cutting Modules feature an excellent cutting quality. The big diameters enable a very time-saving, efficient morcellation.

12,5 mm

15 mm

20 mm

Simple handling –

Change of Cutting Modules / diameters during the operation through the patented click system „Quick Connect“.

On- / Off- Button

For simple handling optionally on the Grip Module or the Foot Pedal.

Variable adjustable –

Controlled speed of Cutting Modules with 100 - 800 rpm constant, powerful efficiency.

Simple processing –

The used materials of the Grip Module and the Trocar Systems enable a standard processing at 134°C.

5 - times - use Check box on the Grip Module to mark the numbers of applications.

Ergonomic Design –


Grip Module lies well in the hand and its weight is balanced.

Compact System –

The small number of components consisting of Operating Unit, Flexible Shaft, Grip Module and accessories are easily connectable to each other and integrate easily in the operation room.

OR-Set Morcellator moresolution

System overview with combination possibilities.



The new moresolution morcellator:
performance, design and functionality perfected.
Superior Lightweight design with More-Power, More-Control
and More-Options.

