



My Tube


Uncomplicated use in gynecology

Reusable and universal instrument to be used in colpotomy, LAVH or TLH.

1 Atraumatic instrument attachment –
In gynecologic-anatomy adapted design

2 Optimal sizes –
Reusable caps in three
sizes

5 Stable access –
One-way valve as an access for
10 mm instruments



3 Simple reprocessing –
robust, biocompatible high performance plastic can be easily processed
at 134 ° C / 3 bar

4 Lightweight –
solid, reusable basic tube

- Good visualization of the uterine vessels for a safe dissection
- Stable pneumoperitoneum during application
- Ease of use, increases patient safety and allows a reduced operating time