

Probe Solution

Vaginal and rectal probes

More security due to a better view

The instrument supports the illustration of the rectal and vaginal walls during the operative repair of the vaginal prolaps.

Ideal for laparoscopic technique: with or without mesh.




- Blunt, atraumatic tips with diameter Ø 28 mm and Ø 40 mm.
- Good visualisation of rectum and vagina for a secure positioning of the seam.
- Usage of high-class materials like polyoxymethylene (POM) and stainless steel.
- Shafts with Ø 8 mm can be used with standard restraint systems.
- Atraumatic tips reduce the trauma.
- Ideal when checking for rectal injury and aids in identification of open blood vessels.
- Easy to handle
- Reusable



Good visualisation of vagina and rectum.


The vaginal/rectal probe serves the pelvic floor reconstruction and has to be inserted vaginally or rectally. For the attachment of mesh implants it allows a better visualisation of vagina and rectum.



Vaginal vault prolapse before operation

Treatment with vaginal/rectal probe and mesh

Atraumatic tips



Ordering information

HDC020	Vaginal probe Ø 40 mm x 35,5 mm, reusable
HDC010	Rectal probe Ø 28 mm x 34,5 mm, reusable